



WE WORK LIFE - *Produtos Farmacêuticos Unipessoal, Lda*
Rua de Campolide Edifício 5 N°351E,8°A1070-034 Lisboa Portugal

CENTRAL WORK - *29 BLVD Tower 1 & 2 - Emaar*
Sheikh Mohammed Bin Rashid BLVD-Downtown Dubai UAE

www.centralworkeau.com

+351 219833081

+351 916655075



Economic Urine drainage bag

- ◆ Medical grade PVC
- ◆ Thickness 0.12 mm
- ◆ Non Return Valve



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB111	1500	Without Outlet	90/130cm & 6.5mm
NLUB112	1500	Cross Valve	90/130cm & 6.5mm
NLUB113	1500	Push Pull Valve	90/130cm & 6.5mm
NLUB114	2000	Without Outlet	90/130cm & 6.5mm
NLUB115	2000	Cross Valve	90/130cm & 6.5mm
NLUB116	2000	Push Pull Valve	90/130cm & 6.5mm

Urine bag with Separate Hanger

- ◆ Medical grade PVC
- ◆ Thickness 0.12 mm
- ◆ Non Return Valve
- ◆ Separate Hanger



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB117	2000	Without Outlet	90/130cm & 6.5mm
NLUB118	2000	Cross Valve	90/130cm & 6.5mm
NLUB119	2000	Push Pull Valve	90/130cm & 6.5mm
NLUB-HANGER	PVC HANGER		

Standard Urine bag with needle free sample port tube clip and bed sheet clamp

- ◆ Medical grade PVC
- ◆ Thickness 0.20 mm
- ◆ Non Return Valve
- ◆ Needle free sample port
- ◆ Tube clamp

02



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB121	2000	Lever Tap	90/130cm & 8/9mm
NLUB122	2000	Cross Valve	90/130cm & 8/9mm
NLUB123	2000	Twist Valve	90/130cm & 8/9mm

Standard Urine bag with PVC Tie

- ◆ Medical grade PVC
- ◆ Thickness 0.18 mm
- ◆ Non Return Valve
- ◆ PVC Tie



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB124	1500/2000	Lever Tap	90/130cm & 8/9mm
NLUB125	1500/2000	Cross Valve	90/130cm & 8/9mm
NLUB126	1500/2000	Twist Valve	90/130cm & 8/9mm

Luxury Urine bag

- ◆ Freeflow drainage
- ◆ Air filter
- ◆ Stronger PVC
- ◆ Bacteria filter
- ◆ Needle-free Sample Port
- ◆ Tube Clamp
- ◆ Bed Sheet Clamp
- ◆ Reinforced Hanger and Drip Chamber

04



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB503	2000/2600/4000ml	Cross/Lever Tap	120cm & 8/9mm

Luxury Urine bag

- ◆ Freeflow drainage
- ◆ Air filter
- ◆ Stronger PVC
- ◆ Bacteria filter
- ◆ Needle-free Sample Port
- ◆ Tube Clamp
- ◆ Bed Sheet Clamp
- ◆ Reinforced Hanger and Drip Chamber



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB504	2000/3000ml	Cross/Lever Tap	120cm & 8/9mm

Luxury Urine bag


- ◆ Freeflow drainage
- ◆ air filter
- ◆ stronger PVC
- ◆ Needle-free Sample Port
- ◆ Tube Clamp
- ◆ Hanger



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB505	2000ml	Cross/Lever Tap	120cm & 8/9mm

Luxury Urine bag

- ◆ Drainage clamp
- ◆ with Holder
- ◆ Tube clamp



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB506	2000ml	Cross/Lever Tap	120cm & 8.8mm

Urine Meter Double Chamber

- ◆ 400 ml urine meter with 2000 ml / 2600 ml capacity collection bag
- ◆ collection bag with needle-free sample port and tube holding clip



06

Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUM101	2000 ml / 2600 ml+400ml	Cross/Lever Tap	120cm & 9/10mm

Urine Meter Single Chamber

- ◆ 500 ml urine meter with 2000 ml / 2600 ml capacity collection bag
- ◆ collection bag with needle-free sample port and tube holding clip



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUM102	2000 ml / 2600 ml+500ml	Cross/Lever Tap	120cm & 9/10mm

Standard Leg Bag with 2 straps



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB611	500/600ml	Pull Push/Cross Valve	10-60cm
NLUB612	750/800ml	Pull Push/Cross Valve	10-60cm
NLUB613	900/1000ml	Pull Push/Cross Valve	10-60cm

Urine Leg Bag with Non woven surface



Item Code	Volume (ml)	Drainage Port	Tube Length and Outside Diameter
NLUB602	500/750ml	Cross/Lever Tap	30-60cm & 8/9mm


Oxygen mask

Product Features:

The prosable oxygen mask without oxygen tube is constructed to supply oxygen or other gases to a patient, and it should be used together with oxygen supplying tube usually. The oxygen mask is made from PVC of medical grade, consists of face mask only

Description

- ◆ Be light in weight, they are more comfortable for patients to wear
- ◆ Universal connector(luer lock) is available
- ◆ Smooth and feathered edge for patient comfort and reducing irritation points
- ◆ The standard length of tube is 7ft, and different length is available
- ◆ CE, ISO approved



Item Code	Description	Tubing
NLOM101	Pediatric Standard/S	with/without Tubing
NLOM102	Pediatric Elongated/M	with/without Tubing
NLOM103	Adult Standard/L	with/without Tubing
NLOM104	Adult Elongated/XL	with/without Tubing


Nebulizer mask

Product Details:

- ◆ Made of Medical grade PVC
- ◆ Adjustable nose clip
- ◆ Elastic strap
- ◆ With 2m tubing
- ◆ Size available:S,M,L,L+,XL
- ◆ Packed in individual PE packing
- ◆ 50pcs/ctn or 100pcs/ctn

Specifications:

- ◆ usually made out of rubber or clear vinyl. It consists of a mask that fits over the patient's nose
- ◆ The standard length of tube is 7ft, and different length is available
- ◆ People with lung conditions often rely on certain medications to help treat these diseases or alleviate symptoms. Oftentimes, a doctor will prescribe a nebulizer for the patient to use, making the medications easier to inhale. In order to receive the medication in this form, a patient normally uses a nebulizer mask



Item Code	Description	Capacity
NLNM101	Pediatric Standard/S	6cc/8cc
NLNM102	Pediatric Elongated/M	6cc/8cc
NLNM103	Adult Standard/L	6cc/8cc
NLNM104	Adult Elongated/XL	6cc/8cc

Nasal oxygen cannula

Product Description:

- ◆ The Nasal Oxygen Cannula is an transporting Oxygen device with double channels, It is used to deliver supplemental oxygen to a patient or person in need of extra oxygen
- ◆ The Nasal Oxygen Cannula is made from PVC in medical grade, consists of connector, mail connected tube, three channel conector, clip, branch connected tube, nostril sucker

Product Features:

- ◆ Green color and white transparent color are both available
- ◆ Transparent clear tube is helpful for inspection of gas delivery
- ◆ Non-toxic PVC is available
- ◆ Softer nose tip can be chosen
- ◆ Available of three basic sizes, L M S, suitable for adult, children and infant use
- ◆ Individual peelable polybag or blister pack Sterile



Item Code	Specification	Qty/ctn	Ctn size(cm)
NLNOC101	Adult Stanrd;Green	250pcs/Ctn	48*40*31
NLNOC102	Pediatric Stanrd;Green	250pcs/Ctn	48*40*31
NLNOC103	Infant Stanrd;Green	250pcs/Ctn	48*40*31
NLNOC104	Adult Stanrd;Trans parent	250pcs/Ctn	48*40*31
NLNOC105	Pediatric Stanrd;Trans parent	250pcs/Ctn	48*40*31
NLNOC106	Infant Stanrd;Trans parent	250pcs/Ctn	48*40*31

Endotracheal tubes

Product Features:

- ◆ Made of non-toxic PVC, the tube is clearer and more smooth
- ◆ High volume and low pressure cuff for choice
- ◆ the main tube is reinforced with stainless steel spring to prevent crushing and kinking
- ◆ With X-ray opaque line and murphy eye

Description

- ◆ An endotracheal tube is a flexible plastic tube that is placed through the mouth into the trachea (windpipe) to help a patient breathe. The endotracheal tube is then connected to a ventilator, which delivers oxygen to the lungs
- ◆ An endotracheal tube may be placed when a patient is unable to breathe on her own, due to a medical emergency, serious illness, or during surgery when a general anesthetic is used
- ◆ Made of medical grade PVC+ superb stainless steel spring
- ◆ With High transparency and lubricated smooth
- ◆ Atraumatic rounded tip(Beveled tip) + Atraumatic murphy eyes(Softly rounded) or Magill type available
- ◆ High volume low pressure cuff and Low profile cuff available
- ◆ Sterile packaging by EO+ X-ray opaque line for accurate placement
- ◆ CE, ISO approved + DEHP FREE available



Item Code	Description	Size	Packing	Outer carton size(cm)
NLET101	Disposable PVC Endotracheal Tube	2.0-10.0mm	1pc/blister, 10pcs/box, 100pcs/ctn	42*38*32cm
NLET102	PVC Reinforced Endotracheal Tube	2.0-10.0mm	1pc/blister, 10pcs/box, 100pcs/ctn	42*38*32cm
NLET103	Silicone Reinforced Endotracheal Tube	5.0-10.0mm	1pc/blister, 10pcs/box, 100pcs/ctn	58*38.5*32.5cm
NLET104	Reusable Silicone Reinforced Endotracheal Tube	5.0-10.0mm	1pc/blister, 10pcs/box, 100pcs/ctn	58*38.5*32.5cm
NLET105	Oral Preformed Endotracheal Tube	2.0-10.0mm	1pc/blister, 10pcs/box, 100pcs/ctn	42*38*32cm
NLET106	Nasal Preformed Endotracheal Tube	2.0-10.0mm	1pc/blister, 10pcs/box, 100pcs/ctn	42*38*32cm

PVC Nelaton Catheter / Suction Catheter/ Stomach Tube / Feeding tube

- ◆ 1. Both frosted and transparent catheters are available
- ◆ 2. Soft, frosted and kink resistant PVC tubing
- ◆ 3. Color coded connector
- ◆ 4. Atraumatic, soft, rounded, closed tip with two lateral eyes for efficient drainage
- ◆ 5. Made from non-toxic, non-irritant, soft medical-grade PVC compatible with Catheter lubricants



Item Code	Description	Specification	Package
NLPC101	Nelaton Catheter	Fr8-Fr24	PE/Blister Pouch
NLPC102	Suction Catheter	Fr8-Fr24	PE/Blister Pouch
NLPC103	Stomach Tube	Fr8-Fr24	PE/Blister Pouch

Silicone Stomach Tube

- ◆ designed for nutrition supplement to stomach and may be recommended for various purposes: for patients who cannot take food or swallow, take enough food by mouth to keep nutrition, congenital defects of mouth, esophagus, or stomach
- ◆ inserted through patient's mouth or nose
- ◆ 1. Be made from 100% silicone
- ◆ 2. Both atraumatic rounded closed tip and opened tip are available
- ◆ 3. Clear depth marks on tubes
- ◆ 4. Color coded connector for identification of size
- ◆ 5. Radio opaque line throughout the tube



Item Code	Description	Specification	Guide Wire
NLSC101	Silicone Stomach Tube	Fr8-Fr24	With or Without

Safety ENFit Enteral Nutrition Feeding Set

Enteral Feeding Bags

- ◆ ① Feeding bag(GravityPump, 500/1000/1200ml)
- ◆ ② PVC tubing(DEHP/Non-DEHP)
- ◆ ③ Roller clamp
- ◆ ④ IV chamber
- ◆ ⑤ Silicon tubing(Optional)
- ◆ ⑥ Silicon connector(Optional)
- ◆ ⑦ Feeding Y connector with cap
- ◆ ⑧ Enfit lock female luer with cap
- ◆ ⑨ Caution tag
- ◆ ⑩ Patient label



ENFit Syringes




Volume	1ml,3ml,5ml, 10ml 20mL,35ml,60ml
Graduation	0.20ml 1.0ml
Applications	Medicine,Flush Hydrate, bolus feed
Optional Package	Bister Pack Back Pack

Suction Connecting Tube With yankauer handle

- ◆ The Connecting tube is usually used together with yankauer handle, and it is intended for suction of the body fluid in combination with aspirator during operation on thoracic cavity or abdominal cavity

Features:

- ◆ 1. Anti-kinking tube to avoid blocking under high pressure
- ◆ 2. Tube length can be determined by clients depending on different requirement
- ◆ 3. EO Sterilization
- ◆ 4. CE, ISO 13485 Certificate approved



Item Code	Description	Type	Air Vent
NLYH101	Yankauer Handle	Round Tip	With or Without
NLYH102	Yankauer Handle	Flat Tip	With or Without
NLYH103	Conneting Tube	180mm	3/16 inch
NLYH104	Conneting Tube	300mm	1/4 inch

Silicone bulb wound drainage reservoir

- ◆ 1. Made of 100% imported medical-grade silicone, with anti-reflux valve
- ◆ 2. Can compatible with any of following drains: Flat Drain, Fluted Drain and Round Drain
- ◆ 3. Squeeze the reservoir and connect with drains for suction, the low-vacuum reservoir can make stable and effective suction for exudates



Item Code	Description	Size
NLWR101	Silicone bulb wound drainage	100ml
NLWR102	Silicone bulb wound drainage	200ml
NLWR103	Silicone bulb wound drainage	400ml

Silicone drainage tube

- ◆ 1. 100% silicone
- ◆ 2. The drainage tube combines a transparent silicone round tube and a white perforated flat drainage tube. The inner wall of the implanted flat segment features tooth-like structure to prevent possible occlusion when the tube is pressed or kinked.
- ◆ 3. 3/4 or full perforation are available for optimal drainage with X-ray opaque effectiveness throughout its entire implanted segment
- ◆ 4. Visual inspection and fluid observation is optimized through the transparent round tube
- ◆ 5. Using trocar for insertion and combining with silicone bulb reservoir for fluids suction and collection
- ◆ 6. Used for various surgeries



Item Code	Description	Type	Size
NLWR104	Silicone drainage tube	Round tube	10-110cm
NLWR105	Silicone drainage tube	Flat tube	25cm
NLWR106	Silicone drainage tube	Cross channel tube	25cm

3 Spring wound drainage

- ◆ 3-springs wound drainage reservoir with 100cm PVC drainage tube and trocar ,bed sheet clamp and robert clamp



Item Code	Description	Type	Size
NLWR111	3 Spring wound drainage	Fr 10 – Fr19	200ml
NLWR112	3 Spring wound drainage	Fr 10 – Fr19	400ml
NLWR113	3 Spring wound drainage	Fr 10 – Fr19	800ml

Accordion wound drainage

- ◆ hollow non-spring wound drainage reservoir with 100cm PVC drainage tube and trocar , bed sheet clamp and robert clamp



Item Code	Description	Type	Size
NLWR114	Accordion wound drainage	Fr 10 – Fr19	400ml
NLWR115	Accordion wound drainage	Fr 10 – Fr19	500ml
NLWR116	Accordion wound drainage	Fr 10 – Fr19	800ml

Medical Grade Self-sealing Pouch

Specification:

- ◆ Brand Name: OEM BRAND
- ◆ Color of film: Blue, Green and White
- ◆ Surface Handling: Flexo Printing
- ◆ Certification: FDA CE, ISO9001, ISO13485, EN868
- ◆ Usage: packaging for hospital, factory and medical device
- ◆ Feature: Disposable
- ◆ Sterilization: Steam and ETO
- ◆ Sealing & Handle: Self Adhesive Seal

Production Description:

- ◆ Medical Grade Paper (60GSM or 70GSM)
- ◆ Transparent Composite Film (White, Blue or Green)
- ◆ Adopt Environmental Printing Ink
- ◆ With STEAM, EO Chemical Indicator
- ◆ Apply to hospital, laboratory, dental clinic nail & beauty supply, tattoo & piercing supply etc



Item Code	Size	Packing
NLSP101	83mm*160mm	200 pcs/box
NLSP102	90mm*260mm	200 pcs/box
NLSP103	135mm*280mm	200 pcs/box
NLSP104	190mm*360mm	200 pcs/box
NLSP105	255mm*390mm	200 pcs/box
NLSP106	200mm*400mm	200 pcs/box
NLSP107	130mm*260mm	200 pcs/box

Medical Grade Pouch Roll

Specification:

- ◆ Brand Name: OEM BRAND
- ◆ Color of film: Blue, Green and White
- ◆ Surface Handling: Flexo Printing
- ◆ Certification: FDA CE, ISO9001, ISO13485, EN868
- ◆ Usage: packaging for hospital, factory and medical device
- ◆ Feature: Disposable
- ◆ Sterilization: Steam and ETO
- ◆ Sealing & Handle: heat-sealing

Product Description:

- ◆ Medical Grade Paper (60GSM or 70GSM)
- ◆ Transparent Composite Film (White, Blue or Green)
- ◆ Adopt Environmental Printing Ink
- ◆ With STEAM, EO Chemical Indicator
- ◆ Apply to hospital, laboratory, dental clinic nail & beauty supply, tattoo & piercing supply etc



Item Code	Size	Packing
NLSR101	50mm*200mm	12 rolls/ctn
NLSR102	100mm*200mm	6 rolls/ctn
NLSR103	250mm*200mm	2 rolls/ctn
NLSR104	400mm*200mm	2 rolls/ctn
NLSR105	75mm*100mm*25mm	8 rolls/ctn
NLSR106	150mm*100mm*55mm	4 rolls/ctn
NLSR107	250mm*100mm*70mm	2 rolls/ctn
NLSR108	400mm*100mm*100mm	2 rolls/ctn