



WE WORK LIFE - *Produtos Farmacêuticos Unipessoal, Lda*
Rua de Campolide Edifício 5 N°351E,8°A1070-034 Lisboa Portugal

CENTRAL WORK - *29 BLVD Tower 1 & 2 - Emaar*
Sheikh Mohammed Bin Rashid BLVD-Downtown Dubai UAE

www.centralworkeau.com

+351 219833081







+351 916655075





Laptop Ultrasound Scanner

- 15 inch LCD screen
- Full digital imaging technology, crystal-clear Image
- Pretty and light
- High Cost-effective
- Broadband multi-frequency probes
- 2 probe connector
- Silica gel backlight gravure laser keyboard
- ARM system powerful function, include report, can work with laser printer
- Built-in battery(can work more than 3 hours)



Arm system, powerful function

Also can customized the color, such as white & gray

Specification:

- Image mode: B,M,2B,B/M,4B
- Display: 15"LCD monitor
- Digital technology: Wide-angle imaging, Panoramic focusing technology, Rolling M technology
- Scanning depth: 300mm
- Imaging Pre-processing: Cloud processing technology, 8-segment TGC, Gain(B&W), Dynamic range, Persistence, M Soften, Noise suppression
- Imaging Post-processing: Image enhancement, Gray Map, Colorize Map, Tissue Specific Imaging, Left/right reverse, Up/down reverse
- Cine-loop: 334 frame cine loop memory
- THI: Tissue Harmonic Imaging
- Store: built-in 8G memory
- Print: by video printer, or by laser printer
- Measurement & Calculations:
 - B mode: Distance, angle, area, volume, trace length, distance ratio
 - M mode: Distance, time, slope, HR
- Software packages: Abdomen, OB, gynecology, small Part
- Transducer connector: 2
- Peripheral ports: 1 BNC video port,2 USB ports,1 networking-port,1 power inlet port
- Power supply: AC 100-240V 50/60HZ <80VA.
- Overall size: 362mm×340mm×60mm
- Gray scale: 256 gray, and pseudo color
- Transducer frequency: 2MHz~11MHz
- Zoom: Pan zoom
- Optimize
- Communication: DICOM 3.0
- Report edit function
- Built-in battery: work 3 hours
- Total weight: 3.5Kg

Standard configuration:

One host, one convex probe, one aluminum trunk



Optional:

Transvaginal, linear, rectal, micro-convex probe;
Biopsy; Trolley; Video printer; Laser printer.

